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FINAL

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Feb 8, 2024

Natupulse® TS 5M**PRD-N°. 30846336**

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Product Specification *FINAL* Effective from Feb 8, 2024

Natupulse® TS 5M **PRD-No. 30846336**

® = Registered trademark of BASF group ™ = Trademark of BASF group DAWF-2024-0121 **Revision 0**

NOT FOR REGULATORY PURPOSES

Test parameter	Requirement	Test method
Enzyme Activity (Mannanase)	≥ 40000 TMU/g	
Loss on Drying	≤ 12.0 g/100g	

The product meets the requirements of this specification.

This document is valid without signature.

Except as expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF, this document and any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law. Unless agreed otherwise in writing, you are responsible for confirming that you have retrieved the most current version of this document from BASF. Any information provided herein can be changed at BASF’s sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy any obligations you may have to undertake your own inspections and evaluations. To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE AND BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, INCLUDING, WITHOUT LIMITATION, ANY LIABILITY FOR ANY DIRECT, CONSEQUENTIAL, SPECIAL, OR PUNITIVE DAMAGES RELATING TO OR ARISING THEREFROM, EXCEPT IN CASES OF (I) DEATH OR PERSONAL INJURY, (II) BASF’S OR ITS AGENTS AND ASSISTANTS WILLFUL MISCONDUCT, FRAUD OR FRAUDULENT MISREPRESENTATION OR (III) ANY MATTER IN RESPECT OF WHICH IT WOULD BE UNLAWFUL FOR BASF TO EXCLUDE OR RESTRICT LIABILITY UNDER THE APPLICABLE LAW.

If you have any further questions or need additional support, please contact your BASF sales representative.

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Method of Analysis FEN/0004/01

Title: Photometric assay of mannanase activity in feed additives

Effective date: April 20th, 2018

Prepared by: Eric Rückert

 20. APR. 2018
Signature, date

Reviewed and approved: Dr. Benjamin Peikert

 20. APR. 2018
Signature, date

Photometric assay of mannanase activity in feed additives

1 Keywords

1.1 To be determined

Mannanase activity

1.2 Matrix

Feed additives such as

- Mannanase SD powder
- Mannanase CC
- Mannanase 8000
- Mannanase 8000 L

2 Principle

The assay procedure is specific for endo-1,4- β -D-mannanase activity. On incubation of dyed carob galactomannan with β -mannanase, the substrate is depolymerised by an endo-mechanism to produce low-molecular weight dyed fragments which remain in solution on addition of ethanol to the reaction mixture. High-molecular weight material is removed by centrifugation, and the colour of the supernatant is measured. By reference to a standard curve, endo- β -mannanase in the assay solution is determined.

3 Working range

The working range of this method is 0.8 - 13.4 TMU per ml.

4 Safety note

The method described involves the handling of hazardous substances. Attention is therefore drawn to the various provisions governing the handling of potentially dangerous material.

5 Reagents

5.1 Azo-Carob-Galactomannan (S-ACGLM), Megazyme

5.2 Citric acid monohydrate p.a., e.g. Riedel-de Haën article 33114

5.3 Sodium hydroxide solution 50 % pure, e.g. Bernd Kraft article 05395.3000

5.4 Tween 20 Fisher BioReagents, e.g. Fisher Scientific article BP337-500

Photometric assay of mannanase activity in feed additives

5.5 Tween 20 solution, $\beta = 100 \text{ g / l}$

Weigh 10.0 g of Tween 20 (5.4) into a 100 ml volumetric flask and dissolve in approximately 80 ml ultrapure water (5.6). Make up to volume with ultrapure water (5.6) and mix. The solution will remain stable for 2 weeks when stored at ambient temperature.

5.6 Ultrapure water

5.7 Citrate buffer, pH 3.50

Weigh 8.93 g of citric acid (5.2) into a 100 ml volumetric flask and dissolve in approximately 80 ml ultrapure water (5.6). Add 2.9 g of sodium hydroxide solution (5.3) and dissolve. While stirring adjust the pH to 3.50 ± 0.05 with sodium hydroxide solution (5.3). Make up to volume with ultrapure water (5.6) and mix. The solution will remain stable for 2 weeks when stored at ambient temperature.

5.8 Citrate buffer, pH 3.50, containing Tween 20

Weigh 89.31 g of citric acid (5.2) into a 1000 ml volumetric flask and dissolve in approximately 800 ml ultrapure water (5.6). Add 29.0 g of sodium hydroxide solution (5.3) and dissolve. Add 1 ml Tween 20 solution (5.5). While stirring adjust the pH to 3.50 ± 0.05 with sodium hydroxide solution (5.3). Make up to volume with ultrapure water (5.6) and mix. The solution will remain stable for 2 weeks when stored at ambient temperature.

5.9 Citric acid solution, $c = 1 \text{ mol / l}$

Weigh 21.0 g of citric acid (5.2) into a 100 ml volumetric flask and dissolve in approximately 80 ml ultrapure water (5.6). Make up to volume with ultrapure water (5.6) and mix.

5.10 Azo-Carob-Galactomannan substrate (2 g / 100 ml)

Carob galactomannan is dyed with Remazolbrilliant Blue R to an extent of approximately one dye molecule per 30 sugar residues and partially depolymerized. The substrate is supplied as powder.

Add 2.0 g of powdered substrate (5.1) to 80 ml of hot ($85 - 90^\circ\text{C}$) ultrapure water (5.6) and stir vigorously on a hot-plate stirrer (6.6). Turn the heat off and continue stirring until the substrate completely dissolves (about 2 h). After the solution has reached room temperature, add 10 ml citrate buffer (5.7) and adjust the pH to 3.50 ± 0.05 with citric acid solution (5.9). Make up to 100.0 ml with ultrapure water (5.6) and mix. Store this solution at 4°C between use. The solution is stable for one week. Shake the solution container before removing aliquots for assays. Because the solution is viscous, it should preferably be dispensed with a dispenser. Before dispensing, the substrate should be warmed to temperature and thoroughly mixed by vigorous shaking.

Photometric assay of mannanase activity in feed additives

5.11 Precipitant solution

Industrial methylated spirits (IMS; 95 % v / v) or ethanol (95 % v / v)

5.12 Thermostable endo-mannanase standard

Use a highly concentrated, water soluble thermostable endo-mannanase preparation with an officially assigned activity. Store the stock in the freezer. Store an amount for daily use in the refrigerator.

5.13 Thermostable endo-mannanase control sample

Use a water soluble thermostable endo-mannanase preparation with an officially assigned activity. Store the stock in the freezer. Store amount for daily use in the refrigerator. If a liquid preparation is used as control sample, store the stock in the refrigerator, too.

6 Apparatus**6.1 Glassware****6.2 Precision balance, with at least 0.01 g readability****6.3 Analytical balance, with at least 0.1 mg readability****6.4 pH meter, with 2 decimal place accuracy****6.5 Bottle-top dispensers, different volumes****6.6 Magnetic hot plate stirrer****6.7 Magnetic stirrer****6.8 Glass fiber filter, e.g. Whatman™ 1823-110 Grade GF/D Glass Fiber Filter without Binder, diameter: 11 cm, Pore Size: 2.7 µm****6.9 Microliter pipette 1000 µl****6.10 Diluter, provided with 0.250 ml and 2.500 ml syringes****6.11 Vortex mixer for test tubes****6.12 Water bath, maintained at 40.0 ± 0.1°C****6.13 Centrifuge, set to at least 3000 g****6.14 Spectrophotometer, provided with 10.00 mm cell and operating at a wavelength of 590 nm**

Photometric assay of mannanase activity in feed additives

7 Preparation of standard, sample and blank solutions

If the period between dilution and incubation will be > 2 h, an ice bath for cooling the diluted solutions is necessary.

7.1 Preparation of the standard solutions

Weigh from the standard (5.12) accurately to within 0.1 mg an amount corresponding to an activity of 67000 TMU in a 50 ml volumetric flask. Dissolve in citrate buffer containing Tween 20 (5.8) by stirring with means of a magnetic stirrer (6.7). Make up to volume with citrate buffer containing Tween 20 (5.8) and mix. Dilute this solution with means of the diluter (6.10) according to the dilution table below in centrifuge tubes. Mix the diluted solutions by using the vortex mixer (6.11). For the last step of the serial dilutions, prepare the dilution in duplicate in centrifuge tubes (standard and blank).

Table 1 Dilution steps of the standard solutions

Standard solution	Activity to be incubated [TMU / ml]	Standard solution [ml]	Citrate buffer (5.8) to be added [ml]	Diluted standard solution [ml]	Citrate buffer (5.8) to be added [ml]	Total dilution factor
A	0.84 ^a	0.030	2.370	0.050	0.950	80 000 ^b
B	3.35 ^a	0.050	1.950	0.100	0.900	20 000 ^b
C	6.70 ^a	0.100	1.900	0.100	0.900	10 000 ^b
D	10.05 ^a	0.150	1.850	0.100	0.900	6 667 ^b
E	13.40 ^a	0.100	0.900	0.100	0.900	5 000 ^b

^a: The exact concentration should be calculated.

^b: The extraction volume was taken into consideration when the dilution factor was calculated.

7.2 Preparation of the control sample solution

Weigh from the control sample (5.13) accurately to within 0.1 mg an amount corresponding to an activity of 71000 TMU in a 50 ml volumetric flask. Dissolve in citrate buffer containing Tween 20 (5.8) by stirring with means of a magnetic stirrer (6.7). Make up to volume with citrate buffer containing Tween 20 (5.8) and mix.

Dilute the solution with citrate buffer containing Tween 20 (5.8) according to the dilution table below in centrifuge tubes by using a diluter (6.10). Mix the diluted solutions by using the vortex mixer (6.11). For the last step of the serial dilutions, prepare the dilution in duplicate in centrifuge tubes (control sample and blank).

Photometric assay of mannanase activity in feed additives

Table 2 Dilution steps of the control sample solution

Control sample solution [ml]	Citrate buffer (5.8) to be added [ml]	Diluted control sample solution [ml]	Citrate buffer (5.8) to be added [ml]	Total dilution factor
0.100	1.900	0.100	0.900	10 000 ^a

^a: The extraction volume was taken into consideration when the dilution factor was calculated.

7.3 Preparation of the sample solutions

7.3.1 Sample containing a carrier completely soluble in water

Transfer, in duplicate, to the nearest 1 mg, the amounts as specified in table 3 into each of two 100 ml volumetric flasks and dissolve in 90 ml of citrate buffer containing Tween 20 (5.8). Make up to volume with citrate buffer containing Tween 20 (5.8) and mix.

Dilute the solutions with citrate buffer containing Tween 20 (5.8) as specified in table 4 a - c. For the last step of the serial dilutions, prepare the dilution in duplicate in centrifuge tubes (sample and blank).

Table 3

Feed additive	Sample weight used [g]
Mannanase SD powder	0.49 - 0.51
Mannanase CC	0.95 - 1.05
Mannanase 8000 L	0.95 - 1.05

Table 4 a

Mannanase SD powder

Dilution step 1 (V 20)		Dilution step 2 (V 10)		Total dilution factor ^c
0.100 ml sample solution (7.3.1)	1.900 ml citrate buffer (5.8)	0.100 ml diluted sample solution	0.900 ml citrate buffer (5.8)	20 000

Method of Analysis FEN/0004/01

Page 7 of 9

Photometric assay of mannanase activity in feed additives

Table 4 b

Mannanase CC

Dilution step 1 (V 16)		Dilution step 2 (V 10)		Total dilution factor ^c
0.100 ml sample solution (7.3.1)	1.500 ml citrate buffer (5.8)	0.100 ml diluted sample solution	0.900 ml citrate buffer (5.8)	16 000

Table 4 c

Mannanase 8000 L

Dilution step 1 (V 14 2/7)		Total dilution factor ^c
0.070 ml sample solution (7.3.1)	0.930 ml citrate buffer (5.8)	142 4/7

7.3.2 Sample containing a carrier not completely soluble in water

Transfer, in duplicate, to the nearest 1 mg, the amounts as specified in table 5, into each of two 250 ml conical flasks. Use a bottle-top dispenser (6.5) to add 200.0 ml of citrate buffer containing Tween 20 (5.8) to each flask. Stir the contents for 30 minutes using a magnetic stirrer (6.7) to prepare a suspension. Filter some of the suspension through a glass fiber filter (6.8) and discard the first few milliliters of the filtrate.

Dilute the filtrates thus obtained with citrate buffer containing Tween 20 (5.8) as specified in table 6. For the last step of the serial dilutions, prepare the dilution in duplicate in centrifuge tubes (sample and blank).

Table 5

Feed additive	Sample weight used [g]
Mannanase 8000	1.95 - 2.05

Table 6

Mannanase 8000

Dilution step 1 (V 14 2/7)		Total dilution factor ^c
0.070 ml sample solution (7.3.2)	0.930 ml citrate buffer (5.8)	2857 1/7

Photometric assay of mannanase activity in feed additives

8 Procedure

8.1 Standard, control sample and sample solutions

Place the centrifuge tubes containing 1.0 ml of the finally diluted solutions into a water bath at $40.0 \pm 0.1^\circ\text{C}$ and start the stop watch. After exactly 5 minutes add with a dispenser (6.5) to each tube 1.0 ml of substrate (5.10) equilibrated at $40.0 \pm 0.1^\circ\text{C}$. After exactly 25 minutes (20 minutes incubation) stop the reaction by adding 5.5 ml of precipitant solution (5.11) with a dispenser (6.5). Mix the solutions with the vortex mixer (6.11). Allow the tubes to equilibrate to room temperature for 10 minutes, then mix the solutions again and centrifuge with at least 3000 g for 15 minutes. Measure the absorbance of the supernatant solution at 590 nm against ultrapure water (5.6).

8.2 Blank solutions

Add to each tube containing the blank solutions 5.5 ml of precipitant solution (5.11) with a dispenser (6.5). Then add 1.0 ml of substrate (5.10). Mix the solutions with the vortex mixer (6.11). After 10 minutes, mix the solutions again. Centrifuge with at least 3000 g for 15 minutes and measure the absorbance of the supernatant solution at 590 nm against ultrapure water (5.6).

At least, one single blank for each matrix is required. The measurement of the blank solutions could be done during the incubation of the samples.

9 Calculation of mannanase activity

The difference of the absorbance (ΔA_{std}) of the standard (t_{20}) minus the absorbance of the blank (t_0) is calculated, respectively. ΔA_{con} (control sample) and ΔA_{sam} (sample) are calculated accordingly.

A quadratic calibration curve showing ΔA_{std} as a function of enzyme activity in TMU / ml is used to determine the enzyme activity.

The equation for the calculation of the mannanase activity of control sample and samples is given below:

$$\text{Mannanase activity} \left[\frac{\text{TMU}}{\text{g}} \right] = \frac{-b + \sqrt{b^2 - 4 \times (c - \Delta A)} \times DF}{2 \times a \times Wt}$$

a, b, c	coefficients of the quadratic equation
ΔA	absorbance difference of the control sample (ΔA_{con}) and samples (ΔA_{sam})
DF	total dilution factor
Wt	sample weight used in gram

Photometric assay of mannanase activity in feed additives

10 Measurement uncertainty

The estimated measurement uncertainty determined as 6-fold measurement of Mannanase 8000 is reflected by the 3-fold standard deviation, which is max. 9 %.

11 Comments

- The method of analysis described is based on the method from Megazyme "ASSAY OF *endo-1,4-β-Mannanase* using AZO-CAROB GALACTOMANNAN"
- The calculated absorbance difference of the sample should be close to the absorbance difference of standard c (~0.48). Otherwise the solution needs to be diluted appropriately.
- The incubation time of 20 min could be reduced to 10 min, with the consequence of an increased coefficient of variance. In this case the dilutions of the standards and samples must be adapted accordingly.
- Total sample volume could be reduced, but the proportion of precipitant solution, enzyme solution and substrate must be constant.
- To increase the stability of the substrate solution, overlay with a few drops of toluene. Stored at 4 °C the solution is stable for at least 12 months if contamination with enzyme is avoided.
- The analytical work was performed by Julia Kojer and Thomas Fischer.